

# VPN CONNECTION TO X-RAY ASSOCIATES

CWeb and Centricity Web refer to the GE Healthcare software programs used to access X-Ray Associates, “Picture Archiving Communication System” (PACS).

**Please Note:** Currently only Windows XP, Vista and 7 are supported. It works best with Internet Explorer 7.0 or 8.0 and security set to medium. There is no support for Mac OS.

Step 1: Establishing **VPN** connection (VPN client is already installed. If not, follow the installation guide)

1. Please click on the **START** button on the bottom left corner and then click on **PROGRAMS** link and find a folder called **CISCO**. Click on the **Cisco Anyconnect VPN client**.




Figure 1: Open the program Cisco Anyconnect VPN client from the Programs menu.

2. A window will pop up like figure 2. Enter **24.244.218.194** in the **Connect to** field and click on **Select**.




Figure 2: Enter 24.244.218.194 in the Connect to field and click on Select.

3. Please accept the certificate by clicking on the **ACCEPT** button.



Figure 3: Accept the certification.

4. After accepting the certificate, the following window will appear. Please enter **xrayweb1** for both the username and password. **Please note, both the username and password is case-sensitive.**



Figure 4: Please enter xrayweb1 for both the Username and Password.

5. Please accept the security alert message once it appears on the screen by clicking on **Yes**.



Figure 5: Please click on Yes to continue.

6. An icon will appear on the bottom right hand corner of the screen where the clock is located. It will look like figure 6. As long as the lock image is visible, the VPN connection is established.




Figure 6: New VPN connection icon.